

[back](#)

MOTomed viva1

1 Operation

Velocity

The operation of the MOTomed viva1 is safe and easy.


Users or therapy/nursing staff can directly set all the functions they need for the daily movement therapy: fast or slow, light or heavy, forwards or backwards. Upon pressing the button "start" the device will start to move the legs slowly and gently. At the same time the user has all functions activated which are required for the training and ensure the users' safety (MovementProtector, SpasmControl with automatic change of direction etc.)

Adjustment possibilities for the training:

- Velocity (number of revolutions) during passive training from 0 to 60 revolutions per minute.
- Resistance levels during active training – from very light up to very heavy (1 to 12 Nm) – can be set in 12 subtle nuances.
- The training duration can be set from 0 to 120 minutes or on continuous operation.

Important for the safety of the user:

Velocity and motor power can be adjusted separately.



[recommend this website to a friend](#)

[top of page](#)