

[back](#)[printversion](#) 

MOTomed gracile

1 Operation

The operation of the MOTomed gracile is safe and easy.

Users or therapy/nursing staff can directly set all the functions they need for the daily movement therapy: fast or slow, light or heavy, forwards or backwards.

Large lettering and pictures on the screen give continuous information about the current training or operating support.

Upon pressing the button "start" the device will start to move the legs slowly and gently. At the same time the user has all functions activated which are required for the training and ensure the users' safety (MovementProtector, SpasmControl with automatic change of direction etc.).

On the opposite graphic, you can see some example screenshots from the MOTomed viva2.

Important safety feature: speed and motor power can be adjusted independent of each other.

[recommend this website to a friend](#)

[top of page](#)

