

[back](#)

MOTOMed viva2


3 SmoothDriveSystem

Ease spasticity, tension and stiffness

MOTOMed spasm therapy differentiates between sudden spasmodic cramps and long term muscular tension. In order to ease these different forms of spasticity, different loosening processes are required. The MOTOMed therefore provides the SmoothDriveSystem for long term muscular tension and the SpasmControl with automatic direction change for sudden spasmodic cramps.

The MOTOMed SmoothDriveSystem: loosens long term muscular tension (spasticity = high muscle tone)

Stiff muscles of arms or legs can be loosened by gently moving the limbs in a flowing motion, slowly releasing muscle tension. Using this knowledge, a particularly gentle drive system has been developed for the MOTOMed viva: the MOTOMed SmoothDriveSystem (employing a belt drive system rather than a gear drive).



[recommend this website to a friend](#)

[top on page](#)